MANUAL PROSEDUR
PEMERIKSAAN *Urine Ekskresi Iodium* (UEI)
LABORATORIUM GAKI
FAKULTAS KEDOKTERAN
UNIVERSITAS DIPONEGORO

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Laboratorium GAKI

Referensi

Makhmudov Amir A, Caldwell Kathleen L. *The Challenge of Iodine Deficiency Disorder A Decade of CDC's Ensuring the Quality of Urinary Iodine Procedures Program*. National Center for Environmental Health. 2011

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
**BAGAN ALIR MANUAL PROSEDUR
Pemeriksaan UEI**

**Ammonium Persulfate Digestion with Spectrophotometric
Detection of the Sandell-Kolthoff reaction**

Lampiran 1.

**URINARY IODINE AMMONIUM PERSULFATE
DIGESTION METHOD**

1. Pipette 250 μl standart in duplo of sol.B
(blank, 12.5, 25, 50, 100, 200),
pipette 250 μl sample and control to each tube

↓

2. Add 1.0 ml Ammonium persulfate to each tube, mix gently

↓

3. Heat all tube for 1 hour at 91 – 95°C in dry bath

↓

4. Cool tubes to room temperature

↓

5. Add 3.5 ml Arsenious acid to each tube,
mix by vortex and stand for \pm 15 mnt

↓

6. Add 400 μl Cerric ammonium persulfate to each tube at 30 second interval, mixing
each with a vortex after addition

↓

7. Read its absorbance at 405 nm exactly 30 mnt after addition of cerric ammonium
persulfate to the first tube

Lampiran 2

URINARY IODINE AMMONIUM PERSULFATE DIGESTION METHOD

Reagents :

1. Ammonium persulfate solution

228.2 gr $(\text{NH}_4)_2\text{S}_2\text{O}_8$ in 1 L H_2O (deionized)

Store in darkness

2. Arsenious acid solution

5 gr As_2O_3 } add H_2O to about 500 ml
25 gr NaCl }

add 60 ml H_2SO_4 p.a

heat gently with stirring to dissolve

cool to room temperature, dilute with H_2O to 1 L
store in the dark bottle

3. Ceric ammonium sulfate solution

24 gr Ceric Ammonium Sulfate + 100 ml H_2SO_4 p.a

add 911.5 ml H_2O

make up 24 hours before use
store in the darknes